# HUMAN CARE

Your Life. Your Way.

## FL 180/FL320

User Manual/Produktbeschreibung/  
Gebruikershandleiding/Bruksanvisning

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---

## **User manual (ENG)**

Product description	3
Safety pre-cautions	4
Overview and Definitions FL180	6
Technical Data FL180	9
Technical Data FL320	11
Assembly FL180	13
Assembly FL320	14
Operation Instructions	16
Examples of Use	19
Maintenance, inspection and test	20
Disassembly	21
Trouble shooting	21
Spare parts	22
Recycling	23
Guidance and manufacturer's declaration	26
CE Label	28

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## Product description

ENG

The intended use of the mobile lift is to transfer an individual from one resting surface to another (such as a bed to a wheelchair). Moving a person in a sling over ANY distance is not recommended. The lift is intended for use in health care, intensive care and rehabilitation


*is a warning triangle used for situations which require extra care and attention.*

### **IMPORTANT!**

*Read the instruction guide for both the patient lift and lifting accessories before use. Lifting and transferring a person always involves a certain level of risk. It is important to completely understand and adhere to the contents of the instruction guides. The equipment should only be used by trained personnel.*

*Please contact Human Care in the event of any uncertainties or questions.*

*The user manual is available at our website for download so that reading programs or enlargement can be used.*

## NOTE!

This instruction guide contains important information for the user of the product. All who use the product should review and completely understand and adhere to the contents of the instruction guide. Remember to keep the instruction guide in a place where it is always available to those using the product.

## Warning!

Certain environments and conditions can limit the correct use of the mobile lifts, including: thresholds, unlevel floor

surfaces, various obstacles, and extra-thick carpets. These environments and conditions can cause the wheels of the mobile lift not to roll as intended, possible imbalance in the mobile lift, and increased exertion by the caregiver. If you are uncertain that your care environment fulfils the requirements for correct use of the mobile lift, please contact your Human Care representative for further advice and assistance.

***Please read and follow the safety precautions listed below. The operation and use of this lift is simple and straightforward. Following these few basic safety precautions will make lifting operations easy and trouble free***

## General pre-cautions

- ALWAYS fit the sling according to the instructions provided.
- ALWAYS check the safe working load of the lift is suitable for the weight of the patient.
- ALWAYS carry out lifting operations according to the instructions in the user manual.
- NEVER disconnect or bypass a control or safety feature because it seems easier to operate the lift.
- DO NOT lift a patient with the castor brakes on. Always let the lift find the correct centre of gravity.
- DO NOT attempt to maneuver the lift by pushing on the mast, boom or patient.
- ALWAYS maneuver the lift with handle provided.
- ALWAYS lower the patient to the lowest comfortable position before transfers.
- DO NOT push a loaded lift at speeds, which exceed a slow walking pace (3Kilometers/ hour 0.8meters/second)
- DO NOT push the lift over uneven or rough ground. Particularly if loaded.
- DO NOT attempt to push/pull a loaded lift over a floor obstruction, which the castors are unable to ride over.
- NEVER force an operation/safety control. All controls are easy to use and do not require excessive force.
- DO NOT park a loaded lift on ANY sloping surface.
- DO NOT use electric lifts in a shower.
- DO NOT charge an electric lift in a bathroom or shower room.

## General pre-cautions

- DO NOT lift a patient unless you are trained and competent to do so.
- THIS lift is for patient lifting. DO NOT use it, or allow it to be used, for any other purpose.
- DO NOT bump the lift down steps, loaded or unloaded.
- DO NOT attempt to negotiate a loaded lift on a slope, which exceeds 1:12 (approximately 5 degrees).
- DO NOT attempt to negotiate a slope without a second helper being present.
- ALWAYS check the lift is not charging before moving as electrical connection maybe damaged.
- the lifting accessory is selected appropriately, in terms of type, size, material and design, with regard to the patient's needs;
- the lifting accessory is correctly and safely applied to the patient in order to prevent injuries;
- the latches are intact; missing or damaged latches must always be replaced;
- the sling's strap loops are correctly connected to the sling bar hooks when the sling straps are stretched up but before the patient is lifted from the underlying surface.

***Never leave a patient unattended during a lifting situation!***

## Before using, make sure that:

- the lift is assembled in accordance with the assembly instructions;
- the lifting accessory is properly attached to the lift;
- the batteries are fully charged
- you have read the instruction guides for the lift and lifting accessories;
- personnel using the lift are informed of the correct operation and use of the lift.
- ALWAYS carry out the DAILY CHECK LIST before using the lift.
- ALWAYS familiarize yourself with the operating control and safety features of a lift before lifting a patient.
- DO NOT use a sling unless it is recommended for use with the lift.

FL180/FL320 has been tested by an accredited testing institute and comply with the requirements specified in the Medical Device Directive 93/42/EEC, for Class 1 products.

FL180 /FL320 comply with the requirements in EN ISO 10535, IEC 60601-1, EN 60601-1-2,

***The lift must not be modified under any circumstance. If you have any questions, please contact Human Care.***


Particular care must be taken when using powerful sources of electromagnetic disturbance, such as e.g. diathermy, so that they are not positioned on or nearby the lift. If you have any questions, please consult Human Care.

## Before lifting, always make sure that:


- ALWAYS plan your lifting operations before commencing.
- the lifting accessories are not damaged;
- the lifting accessory is correctly attached to the lift;
- the lifting accessory hangs vertically and can move freely;


***The lift may not be used in areas where flammable mixtures may occur, for example, in areas where flammable goods are stored.***



Lift in storage position

Lift in storage position

Max load (kg)	180 / 397 lbs
L, min-max (cm)	43-173cm/16,9-68,1"
A, min-max (cm)	145-199 /57,1 - 78,3"
B (cm)	114cm/44,9"
B'(cm)	77cm/30,3"
C(cm)	63-111cm/24,8-43,7"
D(cm)	52-99cm/20,5-39"
E(cm)	112cm/44,1"
F(cm)	11,5cm/4,5"
F'(cm)	4,5cm/1,8"
G(cm)	86cm/33,9"
N.W (cm)	35cm/13,8"
Turning diameter	120cm/47,2"
Material	Aluminium
Total Weight	32kg/70,5lbs
Mast	11,5kg/25,3lbs

Legs	20,5kg/45,2lbs
Wheels	Front: 75 mm (3 inch.) twin wheels Rear: 75 mm (3 inch.) twin wheels with brakes
Emergency lowering	Mechanical and electrical
Lifting interval	1320mm/52"
Sound level	<65 dB
Protection class	IPX4
Electrical data:	24 V
Duty cycle:	10% Max 2 min/18 min
Batteries:	lead-acid type, 2.9 Ah. 24 V.
Battery charger:	100-240 V AC 50-60 Hz, max 650 mA.
Lift motor:	Permanent magnet motor with mechanical safety mechanism 24 V, 7.7 A
Motor for basewidth adjustment:	Permanent magnet motor 24 V, 2.2 A
Operating forces of controls: Hand control: 4.2 N Control panel: 5 N	


This device is built for indoor operation.

Type B according to the degree of protection against electric shock.

Class II equipment.

Max load (kg)	320kg/700lbs
L, min-max (cm)	54/59/64-184/189/194cm / 21,3/23,2/25,2"-72,4/74,4/76,4"
A, min-max (cm)	146/151/156-205,5/210,5/215,5cm / 80,9/82,9/84,8"
B (cm)	150cm/59,1"
B'(cm)	92,5cm/36,4"
C(cm)	76,5-119cm / 30,1-46,9"
D(cm)	64-106cm / 25,2-41,7"
E(cm)	130cm / 51,2"
E'(cm)	110/43,3"
F(cm)	15,5cm / 6"
F'(cm)	8,5cm / 3,3"
G1(cm)	54cm / 21,3"
G2 (cm)	73cm / 28,7"
G3 (cm)	47cm / 18,5"
N.W (cm)	61cm / 24,0"
Turning diameter	167cm / 65,7"
Material	Aluminium
Total Weight	50kg / 110lbs
Mast weight	18kg / 39,6lbs

Weight of legs	32kg / 70,5lbs
Wheels	Front: 100 mm(4")twin wheelsRear: 125mm(5") twin wheels with brakes
Emergency lowering	Mechanical and electrical
Lifting interval	1320mm / 52"
Sound level	<65 dB
Protection class	IPX4
Electrical data:	24 V
Duty cycle:	10% Max 2 min/18 min
Batteries:	lead-acid type, 2.9 Ah. 24 V.
Battery charger:	100-240 V AC 50-60 Hz, max 650 mA.
Lift motor:	Permanent magnet motor with mechanical safety mechanism 24 V, 7.7 A
Motor for basewidth adjustment:	Permanent magnet motor 24 V, 2.2 A
Operating forces of controls: Hand control:4.2 N Control panel: 5 N	

This device is built for indoor operation.

Type B according to the degree of protection against electric shock.

Class II equipment.

The packing consists of a strong, custom made box to ensure the safe arrival of the lift. A number of documents are supplied packed with each lift should be kept safely for future reference.

This lift is suitable for patients in the SITTING, SITTING/RECUMBENT and RECUMBENT positions.

## SLINGS

Human Care has a wide range of slings that meets the various requirements and needs that the user needs.

Human Care recommends that only slings from Human Care is used together with the floor lift.

The choice of sling is made in accordance with Human Care testing protocol. The testing must be done individually for each patient and, in accordance with his or her needs.

Consult the user manual of the sling to ensure that the product is used correctly. Please go to: [www.humancaregroup.com](http://www.humancaregroup.com) for more information about the models and sizes offered.

Example of recommended sling models:

Description	Part. no
Basic sling soft	25050 (Junior S-XL)
Basic sling net	25045 (JuniorS-XL)
Sit sling	25050 (JuniorS-XL)
Bath sling	25055 (JuniorS-XL)
Full body sling	25005 (S-XL)
Full body sling high	25010 (S-XXL)
Comfort sling	25060 (JuniorM-XL)
Toileting sling	25080 (JuniorS-XL)
Toileting sling high	25055 (JuniorS-XL)
Higback plusline	25120 (S-XL)

## CARTON CONTENTS


Place the carton in a clear working area and open carefully. The carton contains:

- Documents package
- Hand control (electric only)
- Battery (electric only)
- Battery charger (electric only)
- Handle component (with screws)
- Sling bar component

**SAFETY NOTE**—the lift is heavy and will need to be lifted with care. You may need assistance to lift the lift from the carton.


## ASSEMBLY INSTRUCTIONS FL180

1. Remove all the parts from the box and place them carefully on the floor.



2. Lock the rear castors and place the mast in the base using push handles.

Hold at handle

- Lock the mast using the knobs and release the braked castors.



**SAFETY NOTE—Avoid trapping fingers. Keep fingers away from the end of the mast when inserting into the mast slot. Tighten the mast-locking device, and unlock the brakes castors at the rear.**

- Release the emergency stop by turning button as indicated on button.
- Before initial use, charge batteries for approx. 4 hours.

## ASSEMBLY OF HANGER BAR

- Take the hanger bar from the box, and test if the clip at the bar ends function well.
- Install the hanger bar at end of the boom.


NOTE—ENSURE the hanger bar is securely installed before usage.

## ALWAYS CHECK THE FOLLOWING BEFORE OPERATION


- The mast is fully locked into position
- The legs of the lift open and close satisfactorily.
- The red emergency stop button, located on the rear of the control box, is in the OFF (out) position.
- Push up and down buttons on the hand control and confirm the boom rises and lowers.
- Fit the battery pack to the lift and make sure the latch holding the pack in place is fully engaged. “Click” in place.

## ASSEMBLY INSTRUCTIONS FL320


11. Remove all the parts from the carton and place on the floor, taking care to protect the finish from damage. Lift the lift carefully from the carton. Stand the lift in the upright position. See below pictures.


12. Assemble the mast and actuator. Follow the above steps to open the lift before installing the actuator.

13. Lay the lift on the floor (ensure the rear castors are securely locked).


4. Drop the mast assembly upright (using the top of the push handle), and drop the mast into the slot, tighten the mast lock knobs, release the braked castors.

**SAFETY NOTE—Avoid trapping fingers. Keep fingers away from the end of the mast when inserting into the mast slot. Tighten the mast-locking device, and unlock the brakes castors at the rear.**

5. Open legs fully. ENSURE leg-locking knob inserts flush against the base casting and the legs will no longer moves freely. The legs will then only operate by the use of foot actuator located at the rear of the base.

6. Release spreader bar from retaining cradle, and the lift is ready for operation.



ENSURE the cradle is securely locked. See picture below.

## ASSEMBLY OF HANGER BAR


1. Take the hanger bar from the box, and test if the clip at the bar ends function well.
2. Install the hanger bar at end of the boom.

NOTE—ENSURE the hanger bar is securely installed before usage.

## ALWAYS CHECK THE FOLLOWING BEFORE OPERATION

- The mast is fully locked into position
- The legs of the lift open and close satisfactorily.
- The red emergency stop button, located on the rear of the control box, is in the OFF (out) position.
- Push up and down buttons on the handset and confirm the boom rises and lowers.
- Fit the power pack to the lift and make sure the latch holding the pack in place is fully engaged. “Click” in place.

## 1. Leg adjustment

The legs on this lift are adjustable for width. The legs can be opened to enable access around armchairs or wheelchairs. For transferring and negotiating narrow doorways and passages the lift legs should be in the closed position.

To achieve the adjustment, push the buttons on the hand control. The adjustment can be carried out with the patient in the lift, but whether loaded or unloaded the adjustment should NOT be made when the lift is moving.


Open the legs

Close the legs

## 2. Castors and braking


The lift has two braked castors, which can be applied for parking. When lifting, the castors should be left free and un-braked, so that the lift will be able to move to the centre of gravity of the lift. DO NOT apply the brakes because if the brakes are applied the patient might swing to the centre of gravity and this may prove disconcerting and uncomfortable.

### 3. Raising and lowering the boom:

The movement of the boom is achieved by a powerful electric actuator, which is controlled by a simple hand control. The hand control has two buttons with directional allows UP and DOWN. The actuator stops automatically at the limit of travel in both directions.


### 4. Emergency Stop:

The Red Emergency Stop Button is located on the rear of the control box and is activated by pressing in. This will cut all power to the lift and only be reset by twisting the button counterclockwise and releasing.

### 5. Mechanical Emergency Lowering:

In the case of a complete electrical failure the electric actuator is fitted with mechanical lowering device (RED BOSS). This will only operate when the lift is under load. The device must be pulled upwards to activate and a slow decent will commence.

### 6. Electrical emergency lowering

In case of a hand control failure, the electrical emergency lowering function on the control

box can be used. On the control box there are holes. The lifting arm can be lowered or lifted via pressing e.g. a pen in the hole. This is a permitted method of lowering/ lifting (you do not get extended functionality by using these when the battery is low).

### 7. Hand control

The hand control has 4 buttons to operate leg spread and lift arm up and down.


The hand control plugs into a socket at the base of the control box. There is a hook on the rear of the hand control, which allows it to be "parked" on the mast or boom when not in use.

### 8. Battery LED indicator

There are three light emitting diodes:

- Green; There is no need to charge the battery
- Yellow; Battery needs charging
- Yellow; Battery needs charging (Incl. acoustic signal when the hand set is activated).


When one continues to use the application at the lowest battery level the diode will blink and continue to give an acoustic signal. At that time it will be possible to drive the lifting arm down only. Moving the lifting arm up as well as leg spread is not possible

## 9. LED indicator Service

Service indicator will blink (yellow diode) every 12 months (norm EN10535 recommendation) or after 8000 cycles (whichever comes first). The first time that the service symbol is shown the control box will provide an acoustic signal. When it is time for service the service indicator will remain illuminated 2 minutes after use.

Afterwards the control box will power down to save the battery. When the service indicator lights up, the system will still be functional and work as normal.

## 10. Resetting of service interval after service has been carried out.

Resetting of service is done by pressing the up and down buttons on the hand control at the same time for 5 seconds. After pressing the buttons for 5 seconds you will receive an audio signal indicating that the timer is reset. The timer will reset the diodes. In the Service Data Toll read out the “Days since last service” and “cycles since last service” will become zero after resetting and start recounting.

## 11. Slings

The selected sling is attached to the hanger bar hooks; each sling is supplied with instructions. Please study the instructions guide for the lifter and sling before using. Once the correct sling has been selected for the patient, attach it to the hanger bar hooks.

## 12. Charging

Connect the charging cable to the lift and mains. The left LED Indicator on the control box (green) is ON when the charger is connected to the mains power. The right LED indicator on the control box (yellow) is ON when the battery is being charged. When the right LED is OFF, the battery is fully charged. Please note that the lift cannot be operated during charging.

The batteries are located in the power pack and are charged through a socket at the base of the control box.

1. Insert the power supply plug into the charge point at the base of the control box. The plug is inserted with a straight push. DO NOT twist the plug in the socket.
2. Plug the charger mains lead/cord into a suitable mains/power outlet
3. Charging is fully automatic. The LED indicator indicates the status of the charging. Note: Even if the charger is left plugged in for extended periods it will not allow the batteries to “overcharge”.
4. To return the lift to operation, switch OFF the mains supply, pull the charger socket out and disengage from the mains supply.

**WARNING!**

**DO NOT use the lift whilst charging is taking place.**

**Be careful not to trip over the charge lead.**

## WARNING

**Check if the lift is not charging before moving as the electrical connection may be damaged.**

**CAUTION: Keep the batteries fully charged. Place the lift on the charge whenever it is not in use. The charger will not allow the batteries to “overcharge”.**

**CAUTION: Never run the batteries completely flat. If the audible warning sounds, complete the lifting operation in hand and place the lift on charge.**

**CAUTION: Never store the lift for long periods without regular charging throughout the storage period.**

**CAUTION: Always make sure the mains power to the charger is switched off before connecting or disconnecting the charger to or from the lift.**

**CAUTION: Never leave the charger plugged into the lift with the mains power off.**

**CAUTION: Never disconnect the charger plug by pulling on the cable.**

**Do not charge an electric lift in a bathroom or shower room.**

## Examples of Use

ENG

Make sure that the sling is fastened securely on all hooks of the suspension so that the person in the sling does not slip or fall out. When moving the patient, stand to the side of the person you are lifting and keep eye contact with the patient to make them feel safe.

**Unbalanced lifting poses a tipping risk!**


**Only push the lift with the handle when moving the lift.**

Always lift securely with the hanger bar vertically above the patient’s centre of gravity. To obtain this consider the placement of the lift. Make sure that the wheels are unlocked before lifting from the bed, chair, toilet or wheelchair.


### To and from bed:

Placement of the lift when lifting from the bed.:

### To and from chair/toilet/wheelchair:

Placement of the lift when lifting from toilet/ chair/wheelchair:

It is recommended that a thorough inspection and test of the lift and its accessories, slings, etc. are carried out on a regular basis.

## ***Service should be carried out once a year.***

The examination and test should be conducted according to the recommendations and procedures below.


Lifetime info lift: 10 years. (when used and serviced as recommended)

Lifetime info battery: 4 years.

## **CHECK LIST**

1. HANGER BAR: Check the hanger bar for freedom of rotation and swing. Check for wear on the central pivot. Check for firm attachment to the boom.
2. BOOM: Check the attachment of the boom to the mast. Make sure there is only minimal side movement of the boom and the boom is free to rotate on the boom mounting.
3. MAST: Check the operation of the mast-locking device. Make sure the mast fully engages into the socket. Check the bottom actuator unit mounting.
4. ACTUATOR: The actuator should require no maintenance other than checking for correct operation and listening for unusual noise.
5. POWE PACK: Check the function of emergency stop button and emergency down/up.
6. BATTERIES: The batteries are housed in the power pack and should not require maintenance other than regular charging as detailed in the charging instructions. Check that the connections remain clean.
7. LEG ADJUSTMENT: Check the legs operate in both full extensions (inward/outward).
8. CASTORS: Check all castors for firm attachment to the legs. Check for free rotation of the castor and the wheels.
9. CLEANING: Clean with ordinary soap and water and/or any hard surface disinfectant. Harsh chemical cleaners or abrasives should be avoided as these may damage the surface finish of the lift. Avoid wetting any of the electric parts.
10. LOAD TEST: The load test should be carried out in accordance with the manufacture's test procedures.
11. BASE AND WHEELS: Ensure base is even and level (all four wheels are on the floor).
12. SLINGS: Check for wear and fraying.
13. LUBRICATION: Oil pivot joints, including mast and boom connections, pedal assembly, and hanger bar joint (only if required).


1. Lower the mast to the lowest possible position, and then insert the hanger bar into retaining cradle. And close legs tightly.

2. Loosen the knob on the mast, lift the mast and boom assembly, and fold carefully toward the legs. Use the push handle to assist this motion.

**NOTE—if you wish to separate the mast and boom completely from the base and legs, undo the mast locking device fully (remove) and lift the mast and boom completely free from the base. Be careful not to lose the mast lock device's, knob and location spindle.**


3. Place all the parts in the box. Pull the restraining strap sharply upwards and extend over the wheels. (This will keep the legs, mast and boom neatly together during transit and storage).

Step 3

**Note – The emergency button should be activated during transport and storage.**

**NOTE—Care should be exercised when folding/unfolding the lift as there is a possible danger of pinching the fingers etc. Please follow the instructions carefully and ask for assistance if you are unsure of the correct procedure.**

Step 4

## Trouble shooting

### 1. Lift does not move up/down when pressing hand control buttons:

- a. Check that charger cable is not connected to lift and mains
- b. Check that emergency stop button has not been pressed
- c. Check battery status
- d. Check that the hand control cable is connected properly to control box

- e. Please contact Human Care

### 2. Lift base width adjustment does not work

- a. Check that charger cable is not connected to lift and mains
- b. Check that emergency stop button has not been pressed
- c. Check battery status

- d. Check that the actuator cable is connected properly to control box
- e. Please contact Human Care


### 3. Charging function does not work

- a. Check that the charger cable is properly connected to the lift and mains
- b. Check that the battery is properly connected
- c. Check that the emergency stop button hasn't been pressed
- d. If you have an external charger, check if battery can be charged using the external charger
- e. Please contact Human Care

### 4. The patient is stuck in a elevated position

- a. Check Battery capacity
- b. Check that emergency stop hasn't been activated
- c. Check that hand control cable is properly connected
- d. Use electrical emergency lowering at the control box by pushing a pen into the emergency lowering hole
- e. Use mechanical emergency lowering at the actuator
- f. Please contact Human care

Spare parts:	Part. no FL180	Part. no FL320
Battery	70975	71234
Hand control	70974	71231
Hanger bar L	70977	70977
Hanger bar M	70978	70978
Hanger bar S	70979	70979
Power cable (EU)	70980	70980
Power cable (UK)	70981	70981
Power cable (US)	70982	70982
Power cable (AUS)	70983	70983
Hanger bar 4-point		71366

**1: PB batteries**

**2: Waste of Electrical  
and Electronic  
Equipment (WEEE)**

*All other parts: Metals*

The batteries from discarded lifts can be returned to Human Care for recycling.

**NOTE: Batteries and electronics containing lead should be recycled according to (WEEE 2002/96 / EC).**

## Guidance and manufacturer's declaration – electromagnetic emissions

The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should assure that it is used in such an environment.

Emission test	Compliance	Electro-magnetic environment - guidance
RF emissions CISPR 11	Group 1	The Human Care lift uses RF energy only for its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Human Care lift is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

## Guidance and manufacturer's declaration – electromagnetic immunity

The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electro-magnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / Burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines n/a. for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode n/a. for common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Human Care lift requires continued operation during power mains interruptions, it is recommended that the lift be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE UT is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electro-magnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80MHz to 2,5GHz</p> <p>10 V/m 800MHz to 2,5GHz</p>	<p>3 Vrms</p> <p>3 V/m</p> <p>10V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Human Care lifts, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = [3,5/3] \sqrt{P}$ $d = [3,5/3] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [7/3] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflected from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Human Care liftsXX is are used exceeds the applicable RF compliance level above, the Human Care liftsXX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Human Care liftsXX.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the Human Care lifts.

The Human Care lifts are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Human Care lifts can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Human Care lifts as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d = 1,2VP	80 MHz to 800 MHz d = 1,2VP	800 MHz to 2.5 GHz d = 0,7VP
0.01	0.12	0.12	0.07
0.1	0,37	0,37	0,22
1	1,16	1,16	0,7
10	3,67	3,67	2,21
100	11.6	11.6	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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