

HUMAN CARE
Your Life. Your Way.

reddot winner 2024

Cosmos
56100, 56105, 56110, 56115, 56120, 56125

User manual / Bruksanvisning / Benutzerhandbuch /
Gebruikershandleiding / Návod k použití

Content

ENGLISH

User Manual	5
Compliances and Standards	6
Safety instructions	6
Specifications	8
Intended use	8
Product description	9
Packaging label	10
Symbols on the packaging label	10
Product label	11
Symbols on the unit and product label	11
Compatibility	12
Emergency functions	13
Attaching and detaching Cosmos	14
Attaching and detaching Cosmos distance strap	15
Hand control	16
Hand control User Interface	17
Battery maintenance	17
Charger	18
Mounting the Cosmos Charging Holder	19
Using the lift	20
Accessories	20
Overview Accessories	20
Data communication	22
Lifting height	22
Care instructions for the lift and accessories	23
Troubleshooting	24
Guidance and manufacturer's declaration	25

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User Manual

Important!

You must read the user manual for your device prior to usage. Keep this booklet and information handy for future use.

“WARNING!” This symbol is used when particular attention is required.

Lifting and moving a person always entails some risk. You should, therefore, read the user manuals for the lift and the lifting accessories carefully. Always ensure that the accessories you use are intended for use with the lift. As a caregiver, you are responsible for the client’s safety, and you must know whether or not the client will be able to cope with the lifting process.

Lifting or lowering a client diagonally or in an angled movement may injure the client or the caregiver.

A diagonal movement can result in a swing of the lift belt, which can cause the client or caregiver to hit other objects in the lifting situation.

Human Care’s lifts are designed only to raise and lower clients in a vertical direction and never diagonally or through angles. The lift must always be positioned directly above the client when lifting. The client must always be positioned above the surface to which they are lowered.

If the product is used in any way that is not recommended by Human Care, then Human Care will not accept responsibility for any injuries or accidents that may occur.

If no other agreement or approval exists, the product owner is responsible for any combinations using accessories or products from suppliers other than Human Care.

Always contact your supplier if you have any questions or anything is unclear.

Anyone using the equipment must have been given adequate training in using the lift and its accessories.

Clients must not be left unattended and should have a trained individual with them when performing a lift.

Human Care’s products are continuously being developed and updated, and we reserve the right to make design changes without prior notice.

Notice to user/patient in case of a serious incident

Any serious incident in relation to the product should be reported to the local contact, who reports to the manufacturer and the competent authority of the country where the user/patient is established.

Compliances and Standards

Human Care is an ISO 13485:2016 certified Swedish medical device company. The Quality Management System is in compliance with US 21 CFR part 820.

The product is CE marked in accordance with EU regulation 2017/745 (MDR), as class I medical

device. The CE mark is on the product.

This product has been tested in accordance with ISO 10535:2021. IEC 60601-1 (Electrical Safety), IEC 60601-1-2 (EMC), and IEC 60601-1-11 (medical equipment for home care).

Safety instructions


Before using the lift for the first time, make sure that:

- The packaging is not damaged, unintentionally opened or exposed to environmental conditions. If so, please contact your Human Care representative.
- Remove the seal from the lift belt after first reading the user manual.
- You have read and understood the user manual for the lift and all of its accessories.
- The lift has been installed in accordance with the installation instructions by certified personnel.
- The lift accessories are intended for use with the lift and are correctly fitted.
- Before use: charge the battery upon arrival. To ensure battery capacity is maintained, the new ceiling lift must be charged for 4 hours immediately upon arrival before it is put in use or stored. This is because the ceiling lift may have been stored for some time and not been charged regularly during that period. The lift should have at least 30 minutes to adjust to a room climate of 20 degrees Celcius.
- Plug the charger into an easily accessible wall socket. The indicator lights on the hand control indicate that the lift is charging. The green indicator light on the charger shows that the charger is powered.


i Important! The lift is designed exclusively for lifting people, and only one at a time.

Before use, always ensure that:


- Personnel who use the equipment have been given adequate training and in the use of the lift and its accessories.
- The lift should not be used by under-aged persons (< 18 years old).
- Keep the lift away from pets or other animals to avoid unexpected operating or damage.
- Keep the product away from lint, moisture, dust, (sun)light, or other heating sources to avoid overheating or creating warm surfaces. The product should be able to ventilate.
- Keep the product away from aggressive liquids like chlorine because it can damage the product and its safety.
- The lift belt is not twisted or worn and can move freely in and out of the lift.
- Avoid strangulation risk with the lift band, hand control cable, or safety stop strap.

- Lift accessories, e.g., sling, distance strap, and rail roller, show no signs of wear or other damage.
- The lift is correctly suspended in the distance strap, rail roller, and rail.
- To operate the lift, the lift belt needs to have tension, by a weight of 1 kg (for example a hanger bar).
- The right types, sizes, materials, and designs of accessories have been selected with respect to the client's safety and needs.
- The lift is marked with its Safe Working Load (SWL). This should not exceed the rail system's max load capacity.
- The lift accessories have been correctly installed. The sling's strap loops and suspension must be checked once the loops are fully extended when lifting but before the client is lifted from the underlying surface.
- It is important to always check that the loops are at the bottom of all hooks (see picture below) and securely placed below the locking spring.


Correct

Wrong

- When transferring the client, stand to the side of the client and make sure that arms and legs are not obstructed. Keep eye contact with the client to help them feel safe.

- This product can be used in a wet room such as a bathroom. However, the main unit is not allowed to be placed underwater or directly under a water jet.

Correct

Correct

Wrong

Human Care provides suitable training in the safe operation of the lift and other accessories. Further information about slings, rail systems, and other accessories can be found at

www.humancaregroup.com.

To ensure that the lift belt will not be damaged, it is important that the hanger bar is balanced and that the sling is suspended correctly. Never hold or pull the lift belt when moving a client. It is also important that the lift belt is tensioned as it is wound in.


Modification of the product is not permitted without the manufacturer's consent!

This product has been tested extensively to ensure safety, quality, and durability. Modification can influence these essential performance factors of the product. Human Care shall be contacted in case there is a specific request for modification. The technical team needs to approve the modification in written form before it is implemented and applied with a client.

Specifications

Lifting capacity, max.:	
56100 / 56105 (IRC)	150 kg/330 lbs
56110 / 56115 (IRC)	220 kg/484 lbs
56120 / 56125 (IRC)	272 kg/600 lbs


Specifications:		
(WxHxL):	283x154x244mm / 11"x6"x9.6"	
Weight of lift:	6 kg/13 lbs	
Batteries:	25.2V DC 4.8Ah (1x), Lithium-Ion	
Safety lowering:	Electronic & Mechanical	
Class:	Internally powered equipment	
Type:	BF	
IP rating lift:	IP25- Protected against solid foreign objects of 12.5 mm diameter and greater. Protected against water jets.	
IP rating hand control:	IP27- Protected against solid foreign objects of 12.5 mm diameter and greater. Protected against the effects of temporary immersion in water.	
Noise level:	< 48 dB(A)	
Lifting speed:	load till 221 kg/486 lbs	2.6 m/min
	>221kg/486 lbs	2.1 m/min
Length of lift belt:	2150 mm/ 84.6"	
Start/stop:	soft	
User interface hand control:	see dedicated chapter	
Operating forces of hand control	4 N	
Charger:	GlobTek GTM96605-G2A1-R3A (100-240V~, 50-60Hz, 1.5A)	
NRTL:	Yes	
Content in box:	Lift, charger, charging holder, quick reference guide, user manual	

This product is designed for repeated use.

Intended use

The Cosmos ceiling lift provides you with the ability to lift and move a client safely and easily, for example to or from a bed, (wheel)chair or toilet or from the floor. The apparatus can also be used for transfers in bed, walking training, and with the help of

some accessories, for balance training and weighing.
The product is suitable for care situations, for example home care, long-term care, rehabilitation care and institutional care.


Product description

1. Ceiling
 2. Distance strap (accessory)
 3. Lift Cosmos
 4. Quick connect system
 5. Hanger bar (accessory)
 6. Sling (accessory)
 7. Handle
- a. 127-220 mm (depending at choice of rail)
 - b. 100-1000 mm (with distance strap)
 - c. 194 mm (lift cover + anchor key on lift, excluding hanger bar)

Cosmos is a light and compact lift developed with the user, caregiver, and technician in mind. Operating the ceiling lift does not require more than one person.

The quick connection system allows for quick mounting and servicing. The handle makes it easy to carry the lift and to hold during mounting.

Due to its innovative design, Cosmos can be installed as a portable lift that can be removed after use. This means that the same lift can be used in several different rooms or wards. Besides, installation directly on the rail roller is possible for permanent use to have the lift always available


The ceiling lift has been designed and tested for use indoors and is classified as an IP25 product.


The product must not be immersed in water.

Immersion to water can cause a short circuit, resulting in the lift stopping to function. This product can be used in a wet room such as a bathroom. However, the main unit is not allowed to be placed underwater or directly under a water jet; see illustrations in the safety instruction chapter.

The climatic conditions should be as follows: ambient temperature from 5 °C to 40 °C, relative humidity of 15% to 90% non-condensing, and atmospheric pressure range from 700 hPa to 1060 hPa.


Stationary

Portable

Packaging label

The packaging label contains information about the product. This label is located on the package.

Symbols on the packaging label

	CE marking
	Manufacturer
	Read the instructions and make sure you understand them before using this product
	Article number/reference number "-XX" = country specific code (AU-Australia, GB-Great Britan, EU-Europe, JP-Japan, NA-North America)
	Serial number

	Medical Device
	Temperature limit Storage/Transport
	Humidity limit Storage/Transport
	Atmospheric pressure limit Storage/Transport
	Unique Product Code
	Unique Device Identification, including (01) EAN code, (11) production date, (21) serial number

Product label

The product label contains information about the lift. This label is located on the cover of the lift.

Symbols on the unit and product label

	CE marking
	Date of manufacture
	Manufacturer
	Max. load
	Read the instructions and make sure you understand them before using this product
	Article number/reference number "-XX" = country specific code (AU-Australia, GB-Great Britan, EU-Europe, JP-Japan, NA-North America)
	Serial number

	QR-code to user manual.
	Medical Device
	Caution
	The lift contains lithium batteries and electronics and it must be recycled appropriately
	Third party certification marking (not applicable for all regions and versions)
	Type-BF product
	Indoor use only
	Unique Product Code
	Unique Device Identification, including (01) EAN code, (11) production date, (21) serial number

Compatibility

Human Care Group slings are designed for use with all Human Care lifts and hanger bars. The products are also compatible with products of other suppliers utilizing the same attachment method for securing slings to a hanger bar and lift. However, due to the wide range of products worldwide, Human Care cannot be responsible for any errors that may result from improper use or application of combinations of slings, hanger bars, and lifts of other suppliers. The use of the combined system is at the prescriber's own risk and liability.

To ensure the safety of the client and caregiver, Human Care recommends the following minimum requirements:

1. When using slings with lifts and hanger bars, made by manufacturers other than Human Care, the combined system needs to be designed with the same sling attachment method for the hanger bar of floor lifts and ceiling lifts, meaning loop-style slings are combined with a loop-style hanger bar, clip style slings are combined with a clip style hanger bar.
2. The involved patient lifts and slings must be CE marked in accordance with Regulation 2017/745 and comply with standard requirements for patient lifts as stated in ISO 10535.
3. The caregiver must follow the manufacturer's instruction manuals and recommendations with respect to the use, care, client, caregiver, cleaning, and inspection of the slings and lifts.
4. Different products on the assembled lift system: i.e., lift, hanger bar, sling, scale, and other accessories, may have different max. loads. It is always the lowest maximum allowable load indicated for the respective products and the system that applies to the entire system. Always check the markings on all individual products in the system.
5. Individual risk assessment, including a practical test of the system and validation with the specific combination used, are required to affirm that the size of the sling is correct for the intended use and compatible to the width and design of the hanger bar. The combined system needs to be documented by the competent assessor.


If there are any questions, please contact your local Human Care office or distributor.

Emergency functions

Emergency stop


In an emergency, the emergency stop can be activated by pulling the red ribbon. The red button will come out, and the lift will stop.

To reset the emergency stop, push the red button back with a finger or tube. If the emergency stop is not reset, the lift will not charge even if plugged in.

Safety lowering

Firstly, pull on the red ribbon, which activates the emergency stop.

Secondly, pull and hold the red ribbon. This will activate the safety lowering.

Always ensure that safety lowering is performed as safely and comfortably as possible for the client.

The safety lowering must only be used in emergency situations.


The belt speed during safety lowering is faster

than the regular speed using the hand control.

As a result, it is important to maintain control of the client while lowering to avoid injury to the client or caregiver.

The safety lowering function is accessible separately (red ribbon) from the normal up/downward movement which is on the hand control.

In the case of a portable installation, the red belt can be cut at a functional length to avoid strangling risk.

Safety catch

The lift is equipped with a safety catch. This prevents involuntary falls. If the lift belt is pulled out too quickly, the safety catch will engage, thus preventing the fall.


Emergency lowering

In the event of power loss or electrical failure in the product, with the client in the hoist, the mechanical emergency lowering is activated, and the client is lowered in a safe and reliable manner. For children or light clients, the caregiver needs to facilitate pulling down the client.


Suppose a lift is to be used by a client living independently. In that case, some form of communication (e.g., alarm system or phone) device shall be installed in the lift area so that the client can summon assistance in an emergency.

Attaching and detaching Cosmos


Stationary set up - Attach

Cosmos can be mounted directly in the rail roller, located in the rail, with the quick connection system. Hold the lift with the handle facing downwards when connecting the Cosmos to the rail roller. Ensure that the anchor key on the lift's top is parallel with the opening of the rail roller. Insert the lift (1), so the anchor key slides into the opening, and turn the lift 90 degrees (2) to ensure that the locking mechanism has been activated (you will hear a click) and the green mark on the anchor key is visible (3). If a lift with in-rail charging is used - connect the cable.

Stationary set up - Detach


If a lift with in-rail charging is used disconnect the in-rail charging cable. Push the lift upwards towards the rail roller until the green mark on the anchor key is not visible anymore, then turn the lift 90 degrees. The lift then unlocks from the rail roller and can be removed.

Portable set up - Attach

Place the lift on a suitable surface, with the handle facing upward. To let the lift belt out, gently pull by holding the anchor key on the belt and simultaneously press the Down button on the hand control.

When connecting to the rail roller, insert and press the anchor key on the lift belt into the locking device of the rail roller (1) and turn the anchor key 90 degrees (2) and release (3), ensure that the locking mechanism has been activated (you will hear a


click) and the green mark on the anchor key is visible (4).

Portable set up - Detach

Lower the lift and place it on a suitable surface with enough lift belt pulled out to be able to turn the anchor key. Hold the anchor key on the lift belt and press it into the


Attaching and detaching Cosmos distance strap

Distance strap: A distance strap can be used between the rail roller and the lift to adjust to the recommended height. For the available distance straps, see the Overview Accessories section.

Attach the distance strap to the rail roller


When connecting to the rail roller, insert the anchor key parallel to the opening of the rail roller, press the anchor key on the distance strap into the locking device of the rail roller (1), and turn the anchor key 90 degrees (2) and release (3), ensure that the locking mechanism has been activated (you will hear a click) and the green mark on the anchor key on the distance strap is visible (4).

To detach the distance strap, follow the same technique but reverse.

Attach the distance strap to the lift in stationary set up

Put the anchor key holder of the distance strap above the anchor key on the ceiling lift and slide the anchor key holder over the anchor key (1), then turn the anchor key holder (2) to secure the connection and let the ceiling lift hang in the distance strap (3).

To detach the distance strap, follow the same technique but reverse.

Attach the distance strap to the lift in portable set up

Insert the anchor key of the lift belt into the anchor key holder of the distance strap from above (1), and then turn it downwards inside the holder (2) to secure the connection and let the ceiling lift hang in the distance strap (3).

To detach the distance strap, follow the same technique but reverse.

The green mark on the anchor key should always be visible between all possible connections to ensure that the set up is secured.

Hand control

The hand control activates the lift when touching any button. The hand control has an Up button (1) and a Down button (4) for raising and lowering the lift. Press the Up button on the hand control to lift the client. Press the Down button on the hand control to lower the client. Pressing the button activates a soft start, and releasing it causes a soft stop. (2) and (3) are for power traverse model. The hanger bar should be connected because its weight is required to activate the up/down movement of the lift. The hand control is also equipped with information about the status of the lift:

- a. service light
- b. warning triangle
- c. battery status

Do not pull the hand control to move the lift along the rail.



The cable and connection of the cable to the ceiling lift can break due to pull forces, especially over time. Move the lift by guiding the client in the sling or using the handle at the lift.

On the short side of the hand control, there is an input for charging the lift and for communicating with the lift (USB-C).

For people with a pacemaker, the hand control should be kept away from their upper body, with a minimum distance of 153 mm/6 inch, because of the magnetic field.

This minimum distance will avoid disruption of the pacemaker. The hand control should be kept away from the persons upper body during the complete lifting process.

Hand control User Interface

The user interface of the hand control has symbols to inform about the status of the product.

On/Off		Press any key to turn on. The lift is on and ready for use. The lift turns off and goes into sleeping mode when not in use. Press any button to reactivate.
Lighted arrows		The green lighted arrows are the available functions.
Battery symbol		Green light, on: Battery level is full-normal. Green light, pulsating: It might be a good idea to charge up the batteries. The battery symbol also pulsates during charging, but arrow symbols will be off. Green light, fast pulsating: The battery level is critically low. The lift can only lower a client. The lift must be charged up immediately to avoid malfunction. No light: After charging till full battery level, the battery symbol disappears (sleeping mode).
Warning Triangle and Down Arrow		Steady yellow light plus a steady green light in the arrow down indicates overload or overheating, and only the lowering function is available.
Wrench		Steady yellow light: The lift's service interval has been reached. Service and annual inspection must be performed.
Warning Triangle and Wrench		Steady yellow light: Malfunction or the lift's service interval has passed. It is no longer safe to use the lift. Service must be performed immediately. Advise to contact a certified technician (see troubleshooting table).

Battery maintenance

During the ceiling lift's lifetime, no battery maintenance is required as long as it is functioning normally. The battery needs to be charged regularly according to the instructions. If

the battery performance is no longer functional, please contact your local Human Care office or distributor for battery replacement.

Charger

The lift is internally powered by lithium-ion batteries, which are rechargeable with a charger. The battery and charger are included with the product. The battery symbol on the hand control lights up when the batteries need recharging. Plug the charger cable into the hand control, and charging will start. When the lift is charging, the battery symbol pulsates.

The lift can only be used with a charger recommended by Human Care, article number 90575, GlobTek GTM96605-G2*.***






- Charge the lift regularly, preferably every night.
- The charger must be disconnected when you use the lift. The lift doesn't work when it is charging.
- If the emergency stop is activated the lift cannot be charged. Check that the emergency stop is not activated when charging.

- The batteries can not be overcharged.
- Charging should take place in a dry environment.

When the batteries reach a critically low battery level, it will no longer be possible to perform more lifts (up or down). However, there is still enough battery power for a safety lowering.

Mounting the Cosmos Charging Holder

Insert the USB-C contact, as shown in pictures 1 and 2, on the backside of the Cosmos Charging Holder. Push the contact all the way in so the clips inside the charging holders are behind the USB-C contact.

Use a finger or a hard object and press the contact all the way down so it locks, as shown in picture 3.

Position the cable as shown in picture 4.

Mount the Cosmos Charging Holder on the wall using two M3 screws. The cable should be towards the floor, and the part with the USB contact towards the ceiling. See picture 5.

Connect the hand control by sliding it upwards to connect to the USB-C contact. See picture 6.

Using the lift

Connect the hanger bar to the anchor key by pushing it into the hole and turn it 90 degrees (see the Cosmos Hanger Bar user manual for more information). Check that the locking is safe by feeling and seeing the green marking.

The hanger bar chosen should be appropriate for the client’s needs and the type of sling used. Place the client in the prescribed sling according to the sling instructions and hang the sling in the hooks on the bar. Check that the sling’s lifting straps are properly anchored in all hooks.

Inform the client what is going to happen and

keep eye contact.

Plan your transfer and prepare the environment before lifting.

Keep client as low as possible to a safe surface.

Adjust transfer distance to condition of the client.

Bring the portable lift or hanger bar in front of the trunk of the client, when (de)attaching the sling connections.

Perform always a safety stop before leaving the surface.

Accessories

If the safe working load for any of the accessories, e.g., slings, hanger bars, and any other accessories, is NOT the same as for the lift, then the lowest max. weight given for any one product is the one that shall apply. Always check the label on the lift and lifting accessories or contact Human Care if you have any questions or if anything is unclear.

Rail rollers - Human Care offers two different rail rollers: the Cosmos Rail Roller, which is used in rails from Human Care, and Cosmos Rail Roller Adapter, which is used with other rails. See the installation manual for more information.

Distance Straps - Gives the ability to have the lift available at lower heights. Available in 9 different lengths, from 20 cm to 100 cm. Other lengths available on request.

Hanger Bars - There are three different types

of hanger bars available; 2- point hanger bar (available in three sizes), 3-point hanger bar, and 4-point hanger bar.

Slings - Human Care’s range of slings is wide enough to meet the user’s various requirements and needs. A fitting protocol must be performed for each client individually and in accordance with their needs. Consult the user manual for the relevant sling to make sure the product is used correctly.

Human Care recommends using only Human Care slings with Human Care lifts. Visit www.humancaregroup.com for further information about which models and sizes are offered.

Scale - For weighing patients in a medical environment.

Overview Accessories

Rail connections	Article nr.	Materials
Rail Roller	56200	Painted steel
Rail Roller Adapter	56201	Painted steel
Rail Roller IRC	56202	Painted steel
Connector hook for other rail rollers	56203	Painted steel
Connector loop long for other rail rollers	56204	Painted steel
Connector loop short for other rail rollers	56206	Painted steel

Distance straps	Article nr.	Materials
Distance strap, 20 cm, 7,9"	94042-200	Fabric/Stainless Steel
Distance strap, 30 cm, 11,8"	94042-300	Fabric/Stainless Steel
Distance strap, 40 cm, 15,7"	94042-400	Fabric/Stainless Steel
Distance strap, 50 cm, 19,7"	94042-500	Fabric/Stainless Steel
Distance strap, 60 cm, 23,6"	94042-600	Fabric/Stainless Steel
Distance strap, 70 cm, 27,6"	94042-700	Fabric/Stainless Steel
Distance strap, 80 cm, 31,5"	94042-800	Fabric/Stainless Steel
Distance strap, 90 cm, 35,4"	94042-900	Fabric/Stainless Steel
Distance strap, 100 cm, 39,4"	94042-1000	Fabric/Stainless Steel

Hanger Bars	Article nr	Materials
Hanger Bar 2-point, 35cm/13.8"	56160	Aluminium/ Glass fibre reinforced polyamid
Hanger Bar 2-point, 45cm/17.7"	56170	Aluminium/ Glass fibre reinforced polyamid
Hanger Bar 2-point, 55cm/21.7"	56180	Aluminium/ Glass fibre reinforced polyamid
Hanger Bar 3-point, 48 cm/18,9"	56190	Aluminium/ Glass fibre reinforced polyamid
Hanger Bar 4-point, L:55xW:37 cm/L:21,7xW:14,6"	56195	Aluminium/ Glass fibre reinforced polyamid
4-point spreader bar Small* 46.5x45x49.5 cm/18.3x17.7x19.4"	92441	Steel
Powered 4-point spreader bar Medium* 46.5x59x49.5 cm/18.3x23.2x19.4"	92442	Steel
Powered 4-point spreader bar Large* 51.5x73x49.5 cm/20.2x28.7x19.4"	92443	Steel

**4-point spreader bar needs to be combined with Connector loop for other rail rollers short, art. nr: 56206.*

Slings	Article numbers starting at...	Materials
Amputee Sling	25140.../25040...	Fabric
Basic Sling	25045.../25054.../25070...	Fabric
Bath Sling	25035...	Fabric
Columbus Sling E-clip	250400...	Fabric
Comfort Sling High Soft	25069...	Fabric
Comfort Sling Soft	25064...	Fabric
Full Body Sling	25010...	Fabric
Full Body Sling E-clip	25005C.../25006C.../25010C.../25011C...	Fabric
Full Body Sling Net	25005...	Fabric
General Purpose sling	NSB-103.../NSB-104.../NSB-105.../NSB-109...	Fabric
High Basic Sling	25032.../25037...	Fabric
Highback Plusline	25120...	Fabric
Limb Sling	NSB-010.../25700...	Fabric
Low Basic Sling	2504.../2507.../	Fabric
Multi Sling	25105....	Fabric
Positioning sling	NSB-4...	Fabric

Silhouette Sling	25015.../25018.../25019...	Fabric
Sit Sling	25030...	Fabric
Toileting Sling	25080.../25085...	Fabric
Toileting Sling E-clip	250850.../250900.../250910...	Fabric
Toileting Sling High	25090...	Fabric

Scale	Article nr.
Scale for Cosmos lift	56205
Cosmos Trolley	70347

Data communication

Cosmos is able to communicate with a standard PC (Win 7/Win 8/Win 10/Win 11). This communication makes it possible to obtain statistics and other information from the lift.

To communicate with the ceiling lift, you need the Human Care Data and Service Application software that can be downloaded from the Human Care partner portal.

Connect the lift to your PC with a standard USB-C cable. Start the Human Care Data and Service Application and follow the instructions. The lift must be switched on for communication to be possible. The lift can not be used while connected to a PC.

The following information can be read from the lift.

- Article number, serial number, and production date
- Max. user weight (SWL)
- Software version
- Service dates
- Periodic usage statistics
- Error log

Connected computer equipment must fulfil the requirements of IEC 60950-1 or IEC 60601-1

The latest software for the lift should always be used. You can find the latest version at the Human Care partner portal.

Lifting height

The lifting height corresponds to the difference between the highest and lowest position the lift can reach. This corresponds to the length of the lift belt. The length of the lift belt is 2150mm (84,6").

If you want to extend the lift belt there are distance straps in various lengths available as extras. The lifting height of the lift is not affected by the length of the distance strap. The total height of the system consists of the rail in which the lift is mounted, the rail roller in the rail, a distance strap if relevant and the lift.

Care instructions for the lift and accessories

Care instructions for the lift

- Plan the charging of the product so that battery capacity is available when lifting is required. For example during the night.
- If the battery is not charged regularly and the lift is not to be used for a longer period, activate the emergency stop, otherwise the battery will be discharged.
- Clean the product before use with the next client.
- Wipe the ceiling lift with a slightly dampened cloth with a water solution using a mild detergent. Or with wipes for chemical disinfection (70% ethanol, 45% isopropanol or similar).
- Disinfection must take place according to the protocol of the facility, in terms of method, duration and frequency with a disinfectant approved by the local care authorities.
- NEVER use strong detergents on the ceiling lift. Approved cleaning agents are isopropanol (45% isopropyl), alcohol with surfactant), ethanol (70% ethanol with denaturalization agent) or comparable. These agents or wipes can be used for cleaning and disinfection.
- NEVER immerse the lift in water, or any other fluids to clean the lift.
- NEVER use sharp objects for cleaning the product or accessories.
- The lift belt must be washed and disinfected as instructed above.
- The lift must always be stored or transported lying horizontally.

Make sure not to damage or remove labels when cleaning.

The essential information on the label will not be readable. Therefore, it is important to take good

care of the labels and to keep them in place.

Annual inspection

The lift must be inspected at least once a year by someone with Human Care qualifications and in accordance with Human Care's directive. Pay particular attention to parts that show wear.

The lift's service indicator will indicate when a service needs to be performed.

Daily inspection

- Always check the parts of the lifting system before use. The check includes a visual inspection of the stitching on both sides of the products. There must be no loose seams or other signs of wear to lift belts, slings, and fabric. If there are, replace the part before use.
- Check the ceiling lift for visible signs of damage.
- Check that the lift is correctly mounted in the distance strap or rail roller.

Monthly inspection

Check the lift belts for damage or wear and tear. Ensure you roll out enough of the lift's belt to allow the entire length to be inspected. Contact a retailer if the belt is damaged to replace it before further use.

Service

Installation, service, and maintenance must only be performed by personnel with Human Care qualifications and in accordance with Human Care's directive. The lift shall not be serviced or maintained while in use with a client. Only original spare parts may be used.

Service agreements

Human Care invites you to sign a service agreement for annual inspection and test loading.

Transport and storage

During transport, or when the lift is not going to be used for some time, the emergency stop button should be pulled out; otherwise, the batteries will self-discharge. To ensure the life time of the battery, it needs to be charged every 28 days. The lift must be transported and stored lying horizontally. The climatic conditions should be as follows: ambient temperature from -25 °C to 70 °C, relative humidity range of 15% to 90%, non-condensing and atmospheric pressure from 700 hPa to 1060 hPa.

Expected service life

The lift is expected to function for 11000 lifting

cycles at maximum lifting capacity or for 10 years with correct service and maintenance.

Battery service life: 600 charging cycles (1 cycle ≈ 150 normal lifts/procedures).

The intensity and loading of the lift define the lifetime of the lift belt. The more symmetrical lifting, the lower the lifting weight, and the lower the frequency of lifts, the more you increase the expected service life of the lift.

Recycling

Once the product has reached the end of its technical life, it must be disposed of in constituent parts according to local authorities. Metal parts must be treated as metal waste for optimum recycling. Plastic parts can be sorted as residual waste.

Disposal of the batteries should be in accordance with the local regulations.

Troubleshooting

Description	Action
- Nothing happens when you press the Up or Down button on the hand control. - No indicator lights up.	- The batteries are fully discharged. Charge the lift immediately. - The emergency stop is activated. Deactivate by pushing the red button.
- Nothing happens when you press the Up or Down button on the hand control. - The service symbol () lights up.	Check the hand control. If the hand control is defective, it needs to be replaced. Contact a service technician.
- Nothing happens when you press the Up button on the hand control. - The warning triangle () lights up.	- The lift motor is overheated. Wait a while and try again. - The lift is overloaded. Lower using the emergency button.
- Nothing happens when you press the Up or Down button on the hand control. - The battery symbol () is pulsating.	Charging in progress. Disconnect the charger.
- Nothing happens when you press the Up button on the hand control. - The Down indicator () lights up.	The lift is in its maximum upper position. Press the Down button on the hand control.
- Nothing happens when you press the Down button on the hand control. - The Up indicator () lights up.	- The lift is in its maximum lower position. Press the Up button on the hand control to raise the lift. - Lift is unloaded. Pull gently on the belt or connect a hanger bar and press the Down button on the hand control to lower it.

Guidance and manufacturer's declaration

Guidance and manufacturer's declaration – electromagnetic emissions

The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should ensure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Human Care lift uses RF energy only for its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Human Care lift is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / Burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines n/a. for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV for line(s) to line(s), ± 2 kV line(s) to Earth	± 1 kV for line(s) to line(s), ± 2 kV line(s) to Earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_r (>95 % dip in U_r) for 0.5 cycle 40 % U_r (60 % dip in U_r) for 5 cycles 70 % U_r (30 % dip in U_r) for 25 cycles <5 % U_r (>95 % dip in U_r) for 5 sec	<5 % U_r (>95 % dip in U_r) for 0.5 cycle 40 % U_r (60 % dip in U_r) for 5 cycles 70 % U_r (30 % dip in U_r) for 25 cycles <5 % U_r (>95 % dip in U_r) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Human Care lift requires continued operation during power mains interruptions, it is recommended that the lift be powered from an uninterruptible power supply or battery.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE U_i is the a.c. mains voltage prior to application of the test level.			

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. Improper operation of the device because electromagnetic interference. Use only original Human Care parts.

Guidance and manufacturer’s declaration – electromagnetic immunity			
The Human Care lift is intended for use in the electromagnetic environment specified below. The customer or the user of the Human Care lift should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V 0,15 MHz- 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	3V 0,15 MHz- 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Human Care lifts, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3,5/3] \sqrt{P}$ $d = [3,5/3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/3] \sqrt{P}$ 800 MHz to 2.5 GHz
Radiated RF EM IEC 61000-4-3	Professional healthcare facility environment: 3 V/m 80 MHz- 2.7 GHz 80% AM at 1 kHz	Home care environment: 10 V/m 80 MHz- 2,7 GHz 80% AM at 1 kHz	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^{a)} should be less than the compliance level in each frequency range. ^{b)}
Proximity fields from RF wireless communications equipment IEC 61000-4-3	see table at page 23	see table at page 23	 Interference may occur in the vicinity of equipment marked with the symbol to the left.

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflected from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Human Care lifts is used exceeds the applicable RF compliance level above, the Human Care lifts should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Human Care lifts.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Human Care lifts.

The Human Care lifts are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Human Care lifts can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Human Care lifts as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d = 1.2VP	80 MHz to 800 MHz d = 1.2VP	800 MHz to 2.5 GHz d = 0.7VP
0.01	0.12	0.12	0.07
0.1	0.37	0.37	0.22
1	1.16	1.16	0.7
10	3.67	3.67	2.21
100	11.6	11.6	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Wireless transmission equipment, mobile telephones, etc. must not be used close to the device as it may affect its function, minimum distance 30 cm/12". Particular caution must be observed when using strong sources of

interference (such as diathermy, etc.) to ensure that no cables or electrical components are placed on or near the device. If in doubt, ask the person in charge of the equipment or the supplier.

Proximity fields of RF wireless communication equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (cm/")	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	30/12	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} , ± 5 kHz deviation, 1 kHz sine	2	30/12	28
710 745 780	704-787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	30/12	9
810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	30/12	28
1720 1845 1970	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	30/12	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	30/12	28
5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	30/12	9

Note:

If necessary to achieve the immunity test level, the distance between the transmitting antenna and the overhead lift should be reduced to 1 meter. The test distance of 1 meter is allowed according to IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier wave is modulated using a square wave signal of 50% in duty cycle.

^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz can be used.

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DANISH

The image shows a detailed product label for the Cosmos 272kg scale. The label is white with black text and features several icons and barcodes. At the top left, the model name 'Cosmos 272kg' is printed in a large, bold font. Below it, a 'REF' box contains the number '56120'. To the right of this, there is a small icon of a factory and the text 'Human Care HC AB Årstaängsvägen 21B 117 60 Stockholm, Sweden'. Further right is a QR code and a circular icon of a person holding a book. Below the 'REF' box, a 'SN' box contains the serial number '123456789012'. To the right of this is the text 'info@humancaregroup.com' and 'www.humancaregroup.com'. Below the 'SN' box, the weight '272kg / 600lbs' is displayed next to an icon of a person and a weight. To the right of this is a crossed-out trash can icon and a 'CE' mark. Below the weight, the date code 'YYYY-MM' is shown next to a 'SE' mark. To the right of this is a 'MD' mark. Below the date code, the label number 'Label 40647C' is printed next to a person walking icon. To the right of this is the technical specification '25,2 V DC 4,8 Ah IP25' and a warning triangle icon. Below the label number, the text 'Hand control: IP27' is printed next to a house icon. At the bottom of the label, there are two barcodes: a large one on the left and a smaller one on the right. Below the large barcode is the number '(01)7331599003312(11)yyymmm(21)123456789012'. Below the small barcode is the number '8 81608 80348 0'. At the very bottom, the text 'Duty cycle: 2 min ON/18 min OFF' and 'Charger: GlobTek GTM96605-G2A1-R3A' is printed. To the right of this is the text 'Conforms to ISO10535 and AAMI ES60601-1' and 'Cert. to CAN/CSA Z10535 and CSA C22.2#60601-1'. At the bottom right is the 'ETL' mark and the text 'Intertek 509340'.

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