Fellow Classic











The Fellow Classic is a high quality rollator ideal for both indoor and outdoor use. Safe and reliable, it is available in five different models - standard, low, high, wide and narrow. The soft braking system is ideal for those users with limited hand mobility.



The Fellow Classic features a unique fork construction that makes it extremely easy to maneuver, especially in smaller spaces. The

construction that makes it extremely easy to maneuver, especially in smaller spaces. The easy-to-use folding system makes storage and transport easy and the solid wheels prevent build up of dirt and rust allowing them to operate and move more smoothly.



Ergonomic Grips

Unique Design

The Fellow Classic's features handles that are ergonomically shaped to distribute the weight put on them evenly. Combined with the soft braking system, this ensures that minimal hand pressure is required to lock the wheels in place. The forward facing handles allow the user to walk in a more upright

position which is easy on both the back and shoulders.

Practical Accessories

The Fellow Classic comes with an assortment of accessories to make day-to-day life easier. From a tray to a cane holder, you can easily find an accessory to suit your needs. The grid easily removes for cleaning and maintenance.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

Accessories/Features



TrayPractical for resting lightweight items.
Maximum load 2kg/4.4lb.
Part no. 08806



Mesh basket Ideal for those great shopping adventures. Part no.: 88051 (low), 88052 (narrow, wide, high), 08805 (standard)



Cane holder Simplifies transport of a cane or umbrella. Part no. 10603



Slow-Down Brake
Prevents the rollator from
moving too fast. Part
no. 10625 (set), 88121
[right], 88122 [left]



Soft Bag Soft carry case ideal for personal items. Part no. 08808



Light kit Illuminate your way while walking. Part no. 91091



IV Stand Part no. 08814



Backrest Backrest for extra comfort while seated. Part no. 10668



Reflector Kit Includes 30 selfadhesive reflectors. Part no. 91090



Security Lock
Flex cord that protects
rollator from theft while
unattended.
Part no. 91092



One-Hand Brake Allows control of brakes with one hand. Must be ordered with rollator. Part no. 10677 (right), 10678 (left)



Curb ClimberHelps raise the rollator onto a curb. Part no. 10661

Technical Specifications

| | Standard | Low | Wide | Narrow | High |
|---------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|---------------------------------|
| | | | | | • |
| Part no. | 88009 | 88109 | 88309 | 88509 | 88609 |
| Handle Height (Min - Max) | 77cm - 95cm / 30.3" - 37.4" | 68cm - 81cm / 26.8" - 37.7" | 77cm - 95cm / 30.3" - 37.4" | 77cm - 95cm / 30.3" - 37.4" | 82cm - 100cm / 32.3" - 39.4" |
| Seat Height | 59cm/23.2" | 51cm/19.9" | 60cm/23.6" | 60cm/23.6" | 64cm / 25.2" |
| Overall Width | 63cm / 24.8" | 62cm/24.4" | 69cm/27.2" | 60cm/23.4" | 63cm/24.8" |
| Overall Length | 69cm / 27.2" | 60cm / 23.6" | 69cm/27.2" | 69cm/27.2" | 69cm/27.2" |
| Width between Handles | 44cm / 17.3" | 44cm / 17.3" | 50cm / 19.7" | 44cm / 17.3" | 44cm / 17.3" |
| Wheel Size | 20cm x 4cm / 7.7" x 1.7" |
| Product Weight | 8kg/17.6lb | 7.7kg/16.9lb | 8.7kg / 19.2lb | 8.1kg / 17.8lb | 8.1kg/17.8lb |
| Maximum User Weight | 150kg/330.7lb | 140kg/308.6lb | 150kg/330.7lb | 130kg/286.6lb | 130kg/286.6lb |