# Staffan RA











The Staffan RA is a highly maneuvrable rollator primarily for indoor use. It is only 59cm (23.2") wide, making it an ideal option for moving around tighter spaces in and around the home.

#### **Compact Design**

The Staffan RA features fixed rear wheels and brakes for better stability and security as well as a narrow overall width that makes it ideal for use in the home. It's easy-to-fold frame makes transport and storage easy. Designed for users that have to deal with walking while suffering from arthritic or neurological conditions, the Staffan RA will allow for maximum stability while performing everyday functions.

#### **Unique Handles**

The Staffan RA features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight. The Staffan RA can also be equipped with a onehand brake to allow a user to control hand functions on either the left or right hand side.

#### **Adjustable Forearm Supports**

The Rebel RA features very unique forearm supports that are adjustable in both height and direction. This ensures that the forearm supports can be positioned in the best possible way for user comfort and stability.

HUMAN CARE SWEDEN (HQ) Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

#### HUMAN CARE CANADA

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com

#### HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

#### HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

## **Accessories/Features**



**Reflector Kit** Includes 30 selfadhesive reflectors. Part no. 91090



Basket Included with rollator in Canada. Ideal for shopping Part no. 40601



**Cane holder** Simplifies transport of a cane or umbrella. Part no. 40603



**Oxygen Holder** Holds oxygen tube firmly in place. Part no. 40615



Security Lock Flex cord that protects rollator from theft while unattended. Part no. 91092



**Light kit** Illuminate your way while walking. Part no. 91091



### One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Part no. 10653 (low), 10654 (high)

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.

## **Technical Specifications**

	Staffan RA Low 59	Staffan RA High 59
Part no.	42221	42231
Handle Height (Min - Max)	86cm - 108cm / 33.8" - 42.5"	107cm - 128cm / 42.1" - 50.4"
Overall Width	59cm/23.2"	59cm/23.2"
Wheel Size	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"
Product Weight	10.3kg/22.7lb	10.6kg/23.4lb
Maximum User Weight	90kg/198.4lb	125kg/275.6lb