


rollator accessories



Identification Tags



Cup Holder



Trays

Our rollator accessories are manufactured to the highest quality standards and are designed to complement and increase the functionality of our rollators. Safe and easy to assemble they give the user the opportunity to customize their rollator for their own requirements and needs.

Security

Optional soft bag kits are available for the Rebel and Carl Oskar series of rollators. These kits include two bags to offer optimal storage for your daily shopping. Our neXus series of rollators comes complete with a zippered soft bag that can remain on the rollator even in the folded position. All bags conveniently close for added security.

Our identification labels feature a slip out card that can be completed to identify your rollator. They come in two vibrant options, making it easy to see your rollator from a distance.

Convenience

A convenient place to rest your drink or magazine, our trays fit quickly and easily on

most of our rollators. We also offer both cane holders and cup holders to help with those outdoor adventures. Our cane holders attach to the rollator frame and can hold either a cane or umbrella for convenient access when needed. Our cup holder holds a standard size bottle.

Accessories for Safety

We offer a selection of slow-down brakes and one-hand brakes for those users in need of additional braking requirements. Our light kit or reflector kit will aid with the visibility of the user when using the rollator in less bright environments.

Note: not all accessories are sold in all countries.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories



Cup Holder

Holds standard size bottle. Fits all rollators. Part no. 90394



Identification Tag

Two options - Green Part no. 90393, Red Part no. 90392. Fits all rollators.



Cane/Umbrella Holder

Simplifies the transport of the cane. Assorted models for each rollator. Part no. 80703 shown [neXus]




Hard Basket

Ideal for shopping. Holds maximum 4.5kg/10lb. Assorted models for each rollator. Part no. 4510 shown [neXus]



Trays

Practical for resting lightweight items. Assorted models for each rollator. Part no. 80702 shown [neXus]




Oxygen Tank Holder

Holds oxygen tank firmly in place. Not available for all models. Part no. 10615 shown [Carl-Oskar/Rebel series].




Curb Climber

Helps raise rollator slightly to rest onto a curb. Standard on neXus 3. Part no. 10612 shown [Carl Oskar/Rebel series]



Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed. Part no. 10620 shown.



Splash Guard

Protects legs and clothes while walking. For Carl-Oskar only. Part no. 10608




Soft Bag

Assorted models for each rollator. Part no. 08808 shown [Fellow Supreme]



One-Hand Brake

Not available for all models. Part no. 10652 shown.



Tested to all relevant standards.

Technical Specifications

Trays			
Part no.		Additional Information	Product Compatibility
80702		Maximum load 2kg/ 4.4lb. Size: 40cm x 31.25cm / 15.7" x 12.3"	neXus 3, neXus 1
10651		For 57cm models	Carl-Oskar, Rebel, Carl-Oskar RA, Rebel RA, Carl-Oskar TW, Rebel TW
10602		For 62cm, 67cm and 72cm models	Carl-Oskar, Rebel, Carl-Oskar RA, Rebel RA, Carl-Oskar TW, Rebel TW
10664			Parade
08806			Fellow Classic
10667			Fellow Supreme

Cane/Umbrella Holder		
Part no.		Product Compatibility
80703		neXus 3, neXus 1
10603		Carl-Oskar, Rebel, Fellow Classic, Fellow Supreme, Carl-Oskar RA and TW, Rebel RA and TW
10675		Parade

This product is manufactured in accordance with all relevant ISO standards.



Learn More About This Product



Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup







youtube.com/channel/UC7fjINh3zTjRx0P1IDf2Da0

Technical Specifications



Hard Baskets

Part no.		Additional Information	Product Compatibility
4510		Maximum load 4.5kg/ 10lb. Size (LxWxD): 43.2 x 29 x 20cm/ 17" x 11.5" x 8"	neXus 3, neXus 1
Assorted		Part no. 40601 Staffan RA and Staffan TW. All others 10601	Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW, Staffan RA and TW
10662			Parade
Assorted		Part no. 08805 (standard), 88051 (low), 88052 (narrow, high, wide)	Fellow Classic
10666			Fellow Supreme




Soft Bags/Bag Kits

Part no.		Additional Information	Product Compatibility
93280		Maximum load 5kg/11lb 5-Pack Part no. 93281 Included with rollator.	neXus 3, neXus 1
08808		Cannot be used on Fellow Classic low	Fellow Series
91095		Maximum load small bag 2.2kg/5lb, large bag 4.5kg/10lb Fits 57cm and 62cm models	Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW
91096		Fits 67cm and 72cm models	Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW

Slow-Down Brake

Part no.		Additional Information	Product Compatibility
4950		Sold in pairs	neXus 3
10620		Sold in pairs	Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW
10625		Optional individual pieces: Right 88121, Left 88122	Fellow Classic, Fellow Supreme






Curb Climber

Part no.		Product Compatibility
4879		neXus 3 Included with rollator
10612		Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW
10661		Integrated in Fellow Supreme, Accessory on Fellow Classic

Backrest

Part no.		Product Compatibility
90391		neXus 3 Included with rollator
10621		Carl-Oskar, Rebel, 57cm, 62cm or 67cm models, Parade
10655		Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW 72cm
10668		Fellow Supreme, Fellow Classic
P48010		neXus oNe Included with rollator

One-Hand Brake

Part no.		Product Compatibility
10652		Carl-Oskar, Rebel
10653		Carl-Oskar RA and TW, Rebel RA and TW, Staffan RA and TW low, Svea TW
10654		Staffan RA and TW high
10677 Right 10678 Left		Fellow Classic
10669 Right 10676 Left		Fellow Supreme

This product is manufactured in accordance with all relevant ISO standards.



Learn More About This Product
















Follow Us on Social Media



humancaregroup.com

Technical Specifications

Assorted Accessories			
Part no.		Additional Information	Product Compatibility
10608		Splash Guard	Carl-Oskar, Carl-Oskar RA and TW
10615		Oxygen Tank Holder	Carl-Oskar, Rebel, Carl-Oskar RA and TW, Rebel RA and TW (Staffan RA and TW Part no. 40605)
08811		Oxygen Tank Holder	Fellow Classic
10671		Oxygen Tank Holder	Fellow Supreme
10679		Newspaper Strap	Fellow Supreme
08814		IV Stand	Fellow Classic
50609-SVEA		Stabilizer Wheel	Svea
91093		Seat Cover	Carl-Oskar, Rebel, Carl-Oskar RA (with seat), Rebel RA (with seat)
22115 Left 22116 Right		HP Handle	Carl-Oskar, Rebel

Assorted Accessories for all rollators		
Part no.		Additional Information
90392		Red Identification Label (sold only in Canada)
90393		Green Identification Label
91092		Security Lock - includes two keys Not sold in North America
91091		Light Kit - battery operated, LED, mounts with magnet or rubber Not sold in North America
91090		Reflector Kit - includes 30 self-adhesive reflectors Not sold in North America
90394		Cup Holder - holds standard size bottle

This product is manufactured in accordance with all relevant ISO standards.



Learn More About This Product



Follow Us on Social Media

