



User information

Bed Raising Block



The bed raising block is buildable and suitable for both leg and base, every unit built 35 mm (1,38 inches).

No loosen blocks are needed.


Suitable for legs up till 85 mm (3,35 inches).

Can be completed with a ring article number: 895-1, to get a whole covered block.

Article number: 895 *Bed Raising Block*
 895-1 *Ring till the Bed Raising Block*

PRODUCT FACTS

- ABS-plastic
- Ring in plexiglass
- Outer dimension 90x90 mm
- Colour, beige
- Weight, 80 gram (2,82 oz)
- Packed 60 pcs/cart.



Human Care HC Sweden AB
Ilabäcksvägen 5
342 32 Alvesta, Sweden
Tel. +46 (0)472-447 70
Fax. +46 (0)472-122 04
www.humancare.se

Page 1 (1)

2008-05-25
Art. no: 99313-ENG Rev A

HUMAN
CARE
Freedom to move!