



You can securely hold on to, and guide, the client when standing up, sitting down and particularly during transfers with the Alpha Stand-up Supports. When being used, the simple handling can easily be combined with a back-friendly work method and elements from nursing concepts. The Alpha Stand-up Supports can be combined with other aids in various ways, depending on the clients needs and functional requirements. The extremely stable carrying handles provide additional safety. The Alpha Looped Belt, PM-3010, can be used optionally for clients with limited body tension. It enables also especially the safe transfer for clients sitting in very narrow

(wheel)chairs where the flexible support belt of the Stand-up Support Magic cannot be applied to the client from behind.



Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Raising Aid / Magic	Polyamide, latex free	red	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-3015	60 x 20	590 g	(01)04250355330154(10)LOT
PM-3012	60 x 20	550 g	(01)04250355330123(10)LOT

Hygiene and maintenance

The product doesn't require any special maintenance. The outer shell can be machine washed at 60 C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

The inner reinforcing components must be removed!

When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.



Made in EU