



The Open Sliding Mat Caremaster is used primarily for repositioning in bed. The inner surface consists of low friction material. The outer surface consists of carbon-fibre reinforced latex-free material which has breathable and liquid-repellent properties. Thanks to the pulling loops, the Open Sliding Mat Caremaster can also be used for heavy patients.

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Caremaster	PES/Carbon, latex free	blue/red	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-2095	175 x 60	485 g	(01)04250355320957(10)LOT
PM-2097	185 x 70	610 g	(01)04250355320971(10)LOT

Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 90°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.



Made in Germany