



Flexible Children's Sliding Board / PM-5020 (w/o)

Flexible Sliding Board with rounded corners / PM-5030

Flexible Sliding Board / PM-5050



The Alpha Flexible Sliding Board has a very smooth surface in order to minimize resistance due to friction. It is used to transfer clients, e. g. from the edge of a bed to a wheelchair, toilet or car seat. The client doesn't have to be able to stand up, but can make the change in position while sitting, either with or without assistance.

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Flexible Sliding Board	Polyethylene, latex free	black	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-5020	58 x 32	1040 g	(01)04250355350206(10)LOT
PM-5030	74 x 32 with rounded edges	1280 g	(010)4250355350305(10)LOT
PM-5050	74 x 32	1320 g	(01)04250355350503(10)LOT

Hygiene and maintenance

The Sliding Boards can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning. When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.



Made in Germany



The distance between free transfer surfaces using the Sliding Board must not be more than 20 cm (10 cm in the case of PM-5020).