

Transfer Aid Merlin Support Belt Merlin



Transfer Aid Merlin (PMM-1100)

Alpha Transfer Air Merlin is a simple and safe aid for assisting people with limited mobility to stand up and change position. It enables carers to transfer clients in an easier, energy-saving and back-friendly manner and it can also mobilise the remaining strength in the client's arms and legs, if applicable, and also give them a feeling of security. The ingenious design in conjunction with the smooth-running rollers allows a particularly safe and easy handling of the Alpha Transfer Aid Merlin. The use of high-quality, rust-free stainless steel also allows maximum stability even in extremely humid environments; for instance - in and out of showers.



Support Belt Merlin (PM-6160)

Alpha Support Belt Merlin assists with the support of a client for sitting and standing, or moving into a wheelchair, onto the toilet or a chair. During a transfer with the Alpha Transfer Aid Merlin it helps the caregiver to control the movement very effectively and obese clients, unstable or nervous clients can be held more securely. Anxious clients are given a sense of security

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Transfer Aid Merlin	stainless steel	grey	200 kg	3 years	5 °C – 40 °C protect against sunlight
Support Belt Merlin	Polyamide, latex free	red	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PMM-1100	53 x 81 x 117	2180	(01)04250355311009(10)LOT
PM-6160	150 x 21	1190	(01)04250355361622(10)LOT

Hygiene and maintenance

The metal frame of the Transfer Aid Merlin is wiped with a normal dust cloth. You can use common liquid household cleaner. The Transfer Aid Merlin can also be treated with an off-the-shelf wipe disinfectant. Verify compatibility first, if applicable according to the disinfectant information. We recommend aldehyde and alcohol free products such as Trionic®, Mikrobac® Tissues, Mikrobac® forte, Bacillo® 30 Foam or Bacillo® 30 Tissues.

PMM-1100        Made in Germany

The product doesn't require any special maintenance. It can be machine washed at 60°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.

PM-6160        Made in EU