

For pain sensitive, heavy or special care clients



The fitted sheet is an accessory item for the turn sheet standard (item no. PM-1582) or turn sheet with hidden handles (item no. PM-2098). The sliding surface of the fitted sheet is made of an extremely low-friction material. The sliding surfaces of the turn sheet and the fitted sheet on top of each other reduce the frictional forces to a maximum; the effect is the same as with a sliding mat.

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Fitted Sheet	Cotton, natural thickening, latex free	Cream	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-2098	200 x 100	1650 g	(01)04250355315847(10)LOT

Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 90°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.



Made in Germany