

Universal  
Many-  
applications  
Supportive



The Alpha Stand-up Aid Quick is a universal aid with many potential applications. It is a safe transfer aid, which allows a transfer of the patient in the bed or from the bed to the wheelchair. It supports the patient in a pleasant and comfortable manner. The Alpha Stand-up Aid Quick provides support on the patient's back. This application invites the client to lean forward in a physiological manner to facilitate a natural motion. This product has 4 handles on both sides in order to be easy-to-use and back-friendly for the caregiver. It can also be used in combination with our Flexible Sliding Board.

## Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Raising Aid Quick	Polyester	grey/red	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-6080	145 x 11	440 g	(01)04250355360809(10)LOT

### Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 60°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.



MD Made in Germany