



The Alpha Bed Ladder acts as a raising aid and supports the client with a physiological pattern of movement when raising the upper part of the body and finally when reaching an upright seating position. The loop at the end makes it possible to attach the ladder to furniture, for instance a - foot end. It has 8 handles in order to perform the raising in small efficient steps.

## Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
<b>Bed Ladder</b>	Polyester	grey/red	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-6090	170 x 8	440 g	(01)04250355360908(10)LOT

### Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 60°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

When necessary, we recommend wipe disinfection with the product Trionic® or comparable products.



Made in Germany