


Safe  
Comfortable  
Support



The Alpha Bed Wedges are designed to actively reposition people requiring nursing care. They provide comfortable support of the client and can be used individually for all requirements. The bed wedges are particularly useful to position the body at night as well as to provide safe support in a sitting position on the bed.



## Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
<b>Bed Wedges</b>	PU-anti slip	blue	N/A	1 year	-5 °C to 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-8010 Set		800 g	(01)04250355380104(10)LOT
PM-8012 2 pcs.	25 x 22 x 13	200 g	(01)04250355380128(10)LOT
PM-8015	50 x 22 x 13	400 g	(01)04250355380159(10)LOT

### Hygiene and maintenance

The Bed Wedges can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning. When necessary, we recommend wipe disinfection with the product Trionic® or comparable products. The use of bleaching agents or alcohol can permanently impair the properties of the Bed Wedges.



Made in Germany